**Gemeinsamer  
Bundesausschuss**

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#### Literatur

1. Reducing mortality from colorectal cancer by screening for fecal occult blood. Minnesota Colon Cancer Control Study. N Engl J Med. 1993; 328(19):1365-71.
2. The risks of screening: data from the Nottingham randomised controlled trial of faecal occult blood screening for colorectal cancer. Gut. 1999; 45(4): 588–592.
3. Repeated screening for colorectal cancer with fecal occult blood test. A prospective randomized study at Funen, Denmark. Scand J Gastroenterol. 1989;24(5):599-606.
4. False negative fecal occult blood tests due to delayed sample return in colorectal cancer screening. Int J Cancer. 2009;125(4):746-50.
5. Analytical Comparison of Three Quantitative Immunochemical Fecal Occult Blood Tests for Colorectal Cancer Screening, Cancer Epidemiol Biomarkers Prev 2011 20; 1492;

6. Difference in performance of fecal immunochemical tests with the same hemoglobin cutoff concentration in a nationwide colorectal cancer screening program, *Gastroenterology*. 2014;147(6):1317-26)
7. Evaluation of quantitative faecal immunochemical tests for Haemoglobin, GMEC, Guilford, 2013
8. Qualitätssicherung in der Routine ist notwendig und gut realisierbar. *Deutsches Ärzteblatt*, Jg. 111; Heft 19, A842 – A843

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#### Referenzen:

- 1) Church TR, Yeazel MW, Jones RM, et al. A randomized trial of direct mailing of fecal occult blood tests to increase colorectal cancer screening. *J Natl Cancer Inst* 2004;96:770-80.
- 2) Segnan N, Senore C, Andreoni B, et al. Randomized trial of different screening strategies for colorectal cancer: patient response and detection rates. *J Natl Cancer Inst* 2005;97:347-57.
- 3) Kapidzic A, Grobbee EJ, Hol L, et al. Attendance and yield over three rounds of population-based fecal immunochemical test screening. *Am J Gastroenterol*. 2014;109:1257-64.
- 4) Homepage of the World Endoscopy Organization: [http://www.worldendo.org/assets/downloads/pdf/activities/fit\\_reports/gmec\\_fit\\_evaluation\\_report.pdf](http://www.worldendo.org/assets/downloads/pdf/activities/fit_reports/gmec_fit_evaluation_report.pdf); accessed March 30, 2015.
- 5) [https://www.deutschepost.de/content/dam/mlm.nf/dpag/images/b/brief\\_postkarte\\_national/broschuere\\_befoerderung\\_von\\_gefaehrlichen\\_stoffen\\_und\\_gege\\_nstaenden\\_teil\\_1a.pdf](https://www.deutschepost.de/content/dam/mlm.nf/dpag/images/b/brief_postkarte_national/broschuere_befoerderung_von_gefaehrlichen_stoffen_und_gege_nstaenden_teil_1a.pdf)

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Brenner, H., Tao, S. Superior diagnostic performance of faecal immunochemical tests for haemoglobin in a head-to-head comparison with guaiac based faecal occult blood test among 2235 participants of screening colonoscopy. *Eur J Cancer*, 2013; 49 (14): 3049-3054

Graser, A., Stieber, P., Nagel, D., Schafer, C., Horst, D., Becker, C.R., Nikolaou, K., Lottes, A., Geisbusch, S., Kramer, H., Wagner, A.C., Diepolder, H., Schirra, J., Roth, H.J., Seidel, D., Goke, B., Reiser, M.F., Kolligs, F.T. Comparison of CT colonography, colonoscopy, sigmoidoscopy and faecal occult blood tests for the detection of advanced adenoma in an average risk population. *Gut*, 2009; 58 (2): 241-248

**Gemeinsamer  
Bundesausschuss**

Hoepffner, N., Shastri, Y.M., Hanisch, E., Rosch, W., Mossner, J., Caspary, W.F., Stein, J. Comparative evaluation of a new bedside faecal occult blood test in a prospective multicentre study. *Aliment Pharmacol Ther*, 2006; 23 (1): 145-154

Hundi, S., Häug, U., Brenner, H. Comparative evaluation of immunochemical fecal occult blood tests for colorectal adenoma detection. *Ann Intern Med*, 2009; 150 162-169

Park, D.I., Ryu, S., Kim, Y.H., Lee, S.H., Lee, C.K., Eun, C.S., Han, D.S. Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening. *Am J Gastroenterol*, 2010; 105 (9): 2017-2025

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**gabControl Fecal Occult Blood**

**M04 K 11**

(informations transfered from original)

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## Table of Contents

1	General Description .....	4
1.1	Background .....	4
1.2	Epidemiology of Colorectal Cancer (CRC) .....	5
2	Diagnostic Methods of Occult Blood .....	5
2.1	Detecting Occult Blood as Screening Test .....	5
2.2	Benefits of iFOBT as Screen Test .....	6
2.3	Confirmatory Test .....	6
3	Product Description: gabControl One Step iFOBT for the Diagnosis of Occult Blood .....	7
3.1	General Description of the Device, Intended Use and Operation of Device; .....	7
3.2	Product Package and Storages .....	7
3.3	Assay Principle of Operation .....	8
4	Comparison and Sensitivity Studies .....	8
4.1	Performance Evaluation Studies .....	8
4.2	Sensitivity .....	9
4.3	Relative Sensitivity and Specificity to the other leading iFOBT .....	10
4.4	Relative Sensitivity/Specificity to Colonoscopy (The Gold-Standard) .....	10
4.5	Absolute gabControl iFOBT Sensitivity to Hemoglobin .....	12
4.6	Specificity .....	12
4.7	Interference Studies .....	15
4.8	Conclusion: .....	16
5	Read Time Studies .....	17
6	Reproducibility Study .....	18
7	Production Details .....	19
7.1	All Materials; .....	19
7.1.1	Materials; .....	19
7.1.2	PVC; .....	19
7.1.3	Major Antibodies; .....	19
7.1.4	Buffer Solution Size 2 ml .....	19
7.1.5	Major Buffers; .....	19
7.2	Production Process: .....	19
7.3	LOT # assignment .....	20
7.4	Flow chart of the manufacturing process .....	21
7.5	Results of Stability studies .....	22
7.5.1	The Accelerating Stability Study; .....	22
7.5.2	Real Time Stability Study; .....	22
8	References .....	24
9	Version History .....	25

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## List of Tables

Table 4-1: Relative Sensitivity/Specificity of the gabControl iFOBT test in comparison with InSure iFOBT test (n = 750).....	10
Table 4-2: Relative Sensitivity/Specificity of the FOBT with Gold-Standard Colonoscopy (n = 260) .....	11
Table 4-3: Absolute Sensitivity of gabControl iFOBT test to hemoglobin .....	12
Table 5-1: Read Time Study .....	17
Table 6-1: Reproducibility Study .....	18
Table 6-2: Reproducibility Study Based on Study Sites .....	18
Table 7-1: Accelerated Stability Study .....	22
Table 7-2: Real Time Stability Study (3 Lots) .....	23

## List of Figures

Figure 1-1: An x-ray study such as a barium enema can also identify the tumor or cancer. ....	4
Figure 7-1: Manufacturing Process Flow Chart .....	21

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## 1 General Description

### 1.1 Background

In 1997, there are approximately 131,200 new cases of colorectal cancer in the US, resulting in 54,900 deaths. Principal risk factors for colorectal cancer include a history of one of the familial polyposis syndromes, familial cancer syndromes, colorectal cancer in first-degree relatives, or a personal history of ulcerative colitis, adenomatous polyps, or endometrial, ovarian, or breast cancer. Colorectal Cancer (CRC) detected at an early stage can be successfully treated with surgery.

Malignancies and, to a lesser extent, polyps, bleed intermittently as the stool moves past them. Performing an annual fecal occult blood test (FOBT) is one of several recommended options for colorectal cancer screening in the average risk population beginning at age 50. Annual or biennial screening with guaiac-based FOBTs (gFOBT) has been shown in large, randomized trials to have a significant and beneficial effect on colorectal cancer incidence and mortality. However, while the specificity of these tests is generally high, sensitivity is poor. Complicated dietary restrictions prior to testing and sampling instructions may limit patient compliance. Newer, immunochemical FOBTs (iFOBT) are reported to have improved performance characteristics compared to guaiac tests, no dietary restrictions,

Figure 1-1: An x-ray study such as a barium enema can also identify the tumor or cancer.

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4/25

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## 1.2 Epidemiology of Colorectal Cancer (CRC)

Colorectal cancer (CRC) is common in the Western world and is rare in developing countries. A sharp increase is seen in Eastern Europe and Japan. ETIOLOGY: Molecular genetics has shown that accumulation of genetic changes is important in the development of colorectal cancer. Mutations in at least four to five genes are required for the formation of a malignant tumor. Environmental mutagenic factors may determine which susceptible individuals grow carcinomas. Environmental risk factors for colorectal cancer are found in a western diet, rich in fat, meat, and animal protein and low in fiber, fruit, and vegetables. The complex interrelations between food components make it difficult to define the precise role of specific food factors.

## 2 Diagnostic Methods of Occult Blood

Screening OB tests are used to look for disease in people who do not have any symptoms. In many cases, these tests can find colorectal cancers (CRC) at an early stage and greatly improve the chances of successful treatment. Screening tests can also help prevent some cancers by allowing doctors to find and remove polyps that might become cancer. There are several tests used for colorectal cancer (CRC).

### 2.1 Detecting Occult Blood as Screening Test.

Fecal Occult Blood Tests are used to detecting Occult Blood in fecal/stool. FOBTs have sufficient sensitivity to detect bleeding that is not visible in the stool.

There are three classes of FOBTs: guaiac-based tests, heme-porphyrin tests, and immunochemical tests.

1) In guaiac-based tests (gFOBT), the pseudoperoxidase activity of hemoglobin turns the guaiac compound blue in the presence of hydrogen peroxide. If certain substances are in the stool the results may be false positive. Patients usually are warned not to eat red meat or certain fruits and vegetables for 72 hours before testing. Use of aspirin or other nonsteroidal anti-inflammatory drugs (NSAIDs) should be avoided for one week before testing. False-negative results may occur in patients taking vitamin C; patients should be advised not to take vitamin C in the week before testing. Patients are encouraged to eat foods high in fiber for one week before testing to cause more rapid stool transit.

2) The heme-porphyrin test, HemoQuant, measures hemoglobin-derived porphyrins, allowing quantitative measurement of hemoglobin in stool. Its use is limited because it has a higher rate of false-positive results, is expensive, and requires laboratory time. Unlike the guaiac-based tests, dietary peroxidases do not cause false-positive

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5/25

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3) Immunochemical FOBTs (iFOBT), detect intact human hemoglobin. Gabmed's one step iFOBT test is based on immunochromatographic principle rather than the chemistry principle, offers a much higher level of sensitivity (50ng/ml). Stool sample is taken with the provided collection buffer, which makes it easier to use than the guaiac-based FOBT (gFOBT), and consumption of red meat or dietary peroxidases does not produce false-positive results.

## 2.2 Benefits of iFOBT as Screen Test

Three case-control studies of iFOBT screening and colorectal cancer incidence or death have been published by the same corresponding author in Japan. Nakajima et al. (2003) found that cases diagnosed with advanced colorectal cancer were significantly less likely than controls to have been screened within the previous 2 or 3 years. Similarly, Saito et al. (1995, 2000) found that cases of death from colorectal cancer were significantly less likely in those screened versus those not screened. Across studies, the risk reductions ranged from approximately 40 to 60%.

## 2.3 Confirmatory Test

When Occult Blood is detected or when the symptoms warrant, a confirmatory test need to be performed.

**Flexible sigmoidoscopy (flex-sig):** A sigmoidoscope is a slender, lighted tube about the thickness of a finger. It is placed into the lower part of the colon through the rectum. This allows the doctor to look at the inside of the rectum and part of the colon for cancer or polyps. Because the tube is only about 2 feet long, the doctor is only able to see about half of the colon.

**Colonoscopy:** A colonoscopy is a longer version of the sigmoidoscope. It allows the doctor to see the entire colon. If a polyp is found, it may be removed it. If anything else looks abnormal, a biopsy might be done. To do this, a small piece of tissue is taken out through the colonoscopy. The tissue is sent to the lab to see if cancer cells are present.

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6/25

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## **3 Product Description: gabControl One Step iFOBT for the Diagnosis of Occult Blood**

### **3.1 General Description of the Device, Intended Use and Operation of Device;**

The gabControl One Step OB (Occult Blood) Test Kit is a simple one step immuno-chromatographic assay for the rapid, qualitative detection of human occult blood in feces. (iFOBT)

The gabControl One Step OB Test Kit (iFOBT) is designed to detect lower levels of fecal occult blood than standard guaiac tests. The basis of the test is an immuno-chromatographic sandwich capture method which yields results that appear more specific to human occult blood and are easier to interpret than those of guaiac-based devices. In addition, unlike guaiac assays, the accuracy of the One Step OB test is not affected by interfering substances and does not depend on the status of the patient at the time the specimen is taken.

The One Step OB Occult Blood Test employs a unique combination of monoclonal and polyclonal antibodies to selectively identify occult blood in feces with a high degree of sensitivity. Elevated levels of human occult blood as low as 50 ng/ml can be detected.

The "Control Line" is used for procedural control. Goat anti-mouse antibody is used for "Control Line" coating. When the gold conjugate reagent flows to the "Control Line" area, the goat anti-mouse antibody will react with the gold conjugate and develop purple red color regardless whether the sample containing any antibody to OB (hemoglobin) or not.

The purpose of this procedural control is to check whether the reagents are working properly or whether a proper amount of sample volume is applied. When the reagents are working properly and a proper amount of sample volume is applied, the "Control Line" becomes purple red within 1-2 minutes or less. If the "Control Line" does not develop color, any test result is not valid. A positive test result indicates the presence of blood in levels of 50 ng/ml or higher in the feces. A negative test result indicates the absence of blood or the presence of blood in levels below 50 ng/ml in the feces.

### **3.2 Product Package and Storages**

OB test kit contains following items to perform the assay;

- iFOBT test device,
- collection device,
- buffer and
- instructions for use.

The test kit should be stored at room temperature. The expiration date was determined under normal storage conditions.

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7/25

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## 3.3 Assay Principle of Operation

The gabControl one step IFOBT test is a sandwich immunoassay test. The test cassette/strip contains a nitrocellulose membrane strip with an immobilized mouse anti-OB (Hemoglobin) - antibody in the test reaction zone (T). An antibody, a goat-anti-mouse antibody is immobilized in the control reaction zone (C) on the nitrocellulose membrane.

90 to 120  $\mu$ l specimen (buffer mixed) sample (containing the respective analytes, hemoglobin) is added onto the gold conjugate pad. The gold conjugate pad contains mouse anti-hemoglobin-antibody which is coupled with Colloidal Gold. The analytes (i.e. hemoglobin) in the collected specimen coupled with the mouse anti-OB-antibody of the Gold Conjugate Pad thus forming an antibody – antigen – Colloidal Gold complex while the liquid is moving along the membrane and transports these complexes together along the membrane by capillary action supported by the Absorb Pad.

When the antibody – antigen – colloidal gold complexes are transported across the membrane and reach the respective immobilized mouse anti-hemoglobin-antibody on the membrane, they are trapped and will form a sandwich complex consisting of: immobilized antibody – antigen – antibody – Colloidal Gold. Only when the applied specimen sample contains a certain concentration of hemoglobin, the formation of this sandwich complex will result in a visible purple color band in the respective test region of the membrane. In case there is not hemoglobin antigen in the specimen sample, the test region of the membrane will remain colorless (no test band visible).

The liquid continues to move to the control area (C) band on the nitrocellulose membrane. There, this conjugate will form a complex with the immobilized anti mouse antibody on the membrane resulting in the formation of a purple colored control (C) test band. This indicates that the test has been performed correctly.

Please note the symbols used in the following pages:

Symbol "+" means that there is a visible purple red color line;  
 Symbol "-" means that there is no visible purple red color line;  
 Symbol "N" means the test result is negative;  
 Symbol "P" means the test result is positive.

## 4 Comparison and Sensitivity Studies

### 4.1 Performance Evaluation Studies

The below described performance studies were conducted at the manufacturer's site, gabmed GmbH or where indicated. The studies are documented in their respective files according to the Quality Assurance Manual at the manufacturer's laboratory.

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8/25

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## 4.2 Sensitivity

The sensitivity testing of the gabControl iFOBT tests were conducted by 1) a comparison study with a commercially available iFOBT test kits (InSure-Enterix) among high risk individuals, 2) a study to determine the sensitivity and specificity relative to colonoscopy (gold-standard) with high risk individuals and individuals with previously negative colonoscopy outcome. 3) Sensitivity test using commercially available controls.

These studies were performed as a multi-center study in following laboratories (July 2002 to June 2003):

1. Central Laboratory of Hospital A (Rainbow Hospital JiaXing)

This hospital has 650 beds with several medical departments. The central laboratory of the hospital has 25 medical laboratory assistants and analyses on a routine basis approximately 2000 patient samples in the area of coagulation, drug monitoring, hematology, endocrinology, infectious disease, pre-natal analyses, microbiology, hormone analyses and histology.

2. Central Laboratory of Hospital B (ShaoSan Number 2 Hospital, Shaosan)

The hospital has 870 beds with several medical departments. The central laboratory of the hospital has 33 medical laboratory assistants and analyses on a routine basis approximately 3000 patient samples in the area of coagulation, drug monitoring, hematology, endocrinology, infectious disease, pre-natal analyses, microbiology, hormone analyses and histology.

3. Central Laboratory of Hospital C (Jenzhow hospital, Zeijian)

The hospital has 500 beds with several medical departments. The central laboratory of the hospital has 17 medical laboratory assistants and analyses on a routine basis approximately 1000 patient samples in the area of coagulation, drug monitoring, hematology, endocrinology, infectious disease, pre-natal analyses, microbiology, hormone analyses and histology.

4. Central Laboratory of Hospital D (People's traditional hospital, Zhijai)

The hospital has 500 beds with several medical departments. The central laboratory of the hospital has 21 medical laboratory assistants and analyses on a routine basis approximately 1000 patient samples in the area of coagulation, drug monitoring, hematology, endocrinology, infectious disease, pre-natal analyses, microbiology, hormone analyses and histology.

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9/25

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### 4.3 Relative Sensitivity and Specificity to the other leading iFOBT

Here, a comparison study is performed with a commercially available iFOBT test kits (InSure-Enterix) among high risk individuals.

A total of 750 patients from the high risk but asymptomatic group were selected. Fecal/stool specimens were collected using the collection devices provided and stored at 2-25 C and tested within 4 days of collection.

The age range of these 750 individuals was from 50 to 81 years, 56% were female and 44% were male. The collected specimens were tested in a commercially available iFOBT test and the gabControl one step iFOBT test. All comparison studies were conducted blindly, i.e. the laboratory technician did not know which samples was OB positive or negative before testing.

**Table 4-1: Relative Sensitivity/Specificity of the gabControl iFOBT test in comparison with InSure iFOBT test (n = 750)**

		InSure iFOBT		
		+	-	
gabControl iFOBT (50ng/ml)	+	60	3	63
	-	2	685	687
		62	688	750

**Conclusion:** The relative sensitivity of the gabControl iFOBT test is 96.8 % (60/62) and the relative specificity is 99.6% (685/688) when compared to the InSure iFOBT. It is important to note that InSure iFOBT is not the gold-standard for colorectal cancer (CRC).

### 4.4 Relative Sensitivity/Specificity to Colonoscopy (The Gold-Standard)

The ultimate goal of iFOBT test is to detect colorectal cancer (CRC), and colonoscopy is the recognized gold-standard for CRC diagnosis. Thus, this study is set up to obtain gabControl iFOBT test's relative sensitivity and specificity to colonoscopy.

Two groups of individuals are formed for this study.

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10/25

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- 1) The first group consisted of 136 very high risk individuals (higher than 50 years of age with at least one case of colorectal cancer (CRC) in a direct blood relative, who never had a colonoscopy performed but are scheduled to have the colonoscopy). The fecal/stool samples for this group are taken prior their colonoscopy procedures.
- 2) The second group consisted of 124 individuals who had colonoscopy procedure performed with a negative outcome within 6 month.

One gFOBT (Hemoccult II) and two iFOBT (gabControl and InSure) are used for testing with the fecal/stool samples collected from the above two groups. The fecal/stool samples are identified in such a way that the technician who performed the FOBT does know the group identity of the samples.

For the purpose of this study, when a colonoscopy find colorectal cancer (CRC) or adenomas > 1cm, the diagnosis is positive, otherwise is considered to be negative. The positive/negative diagnosis by FOBTs are determined by following product inserts respectively.

**Table 4-2: Relative Sensitivity/Specificity of the FOBT with Gold-Standard Colonoscopy (n = 260)**

		Colonoscopy					
		+			-		
		gabControl iFOBT	InSure iFOBT	Hemoccult-II, gFOBT	gabControl iFOBT	InSure iFOBT	Hemoccult-II, gFOBT
gabControl - iFOBT/InSure-iFOBT/Hemoccult II	+	26	24	18	2	1	7
	-	12	14	20	220	221	215
		38	38	38	222	222	222
		68%	63%	47%	99.1%	99.5%	96.8%
		Sensitivity			Specificity		

**Conclusion:** When compare to Colonoscopy, the Gold-Standard of colorectal cancer diagnosis, gabmed's iFOBT test has a sensitivity of 68%, InSure iFOBT has a sensitivity of 63% and the Hemoccult-II gFOBT has a sensitivity of 47%. While the specificity of gabControl iFOBT is of 99.1%, InSure iFOBT is of 99.5% and Hemoccult-IIgFOBT of 96.8% relative to Colonoscopy.

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11/25

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## 4.5 Absolute gabControl iFOBT Sensitivity to Hemoglobin

A set of 100 OB collection buffer were spiked with a fixed level of hemoglobin such as 15 ng/ml, 25ng/ml, 50 ng/ml and 75ng/ml (20 buffers per each hemoglobin spiked level). The remaining twenty samples were spiked with zero hemoglobin.

**Table 4-3: Absolute Sensitivity of gabControl iFOBT test to hemoglobin**

Concentration of Hemoglobin (ng/ml)	gabControl iFOBT  (visible test and control band on test cassette strip = positive)  (no visible test band, but visible control band on test cassette strip = negative)
0 ng/ml	-
15 ng/ml	-
25 ng/ml	-
50 ng/ml	+
75 ng/ml	+

It is therefore concluded that gabControl iFOBT has an absolute sensitivity to hemoglobin of 50ng/ml, and is one of the iFOBTs with the highest sensitivities

## 4.6 Specificity

The specificity is defined as the probability to have a negative result in the absence of the particular condition. The specificity was assessed by studying specimen samples from OB (hemoglobin) negative subjects containing no other closely related substance. Specimen samples that were negative to OB were spiked with other substance. These samples were tested using the gabControl iFOBT test by a replicate of 10. A sample was classified negative, when no purple color band was visible for the Test Line but the purple color "C" control line being visible within 5 minutes. A sample was classified positive, when both the control and test line were visible within 10 minutes (results see Table 5.3).

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12/25

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Table 5.3. Specificity testing of the gabControl iFOBT test using negative and positive control samples spiked with various substances.

Number of Tests	Samples spiked with following substances	Negative Sample to OB, spiked with	Positive Sample to OB, spiked with
		<b>gabControl iFOBT</b> (visible test and control band on test cassette strip = positive)  (no visible test band, but visible control band on test cassette strip = negative)	<b>gabControl iFOBT</b> (visible test and control band on test cassette strip = positive)  (no visible test band, but visible control band on test cassette strip = negative)
10	Acinetobacter calcoaceticus	All (-)	All (+)
10	Proteus vulgaris	All (-)	All (+)
10	Salmonella typhi	All (-)	All (+)
10	Acinetobacter spp	All (-)	All (+)
10	Staphylococcus aureus	All (-)	All (+)
10	Candida albicans	All (-)	All (+)
10	Neisseria gonorrhoeae	All (-)	All (+)
10	Escherichia coli	All (-)	All (+)
10	Neisseria catarrhalis	All (-)	All (+)
10	Gardnerella vaginalis	All (-)	All (+)
10	Neisseria meningitidis	All (-)	All (+)

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14/25

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## 4.7 Interference Studies

Potentially interfering chemical substances such as pain medication, lipids, bilirubin and glucose were supplemented to clinically defined negative normal specimens and clinically defined positive specimens. These samples were tested using the gabControl iFOBT test by a replicate of 10. A sample was classified negative, when no purple color band was visible for the test line but the purple color "C" control line being visible within 10 minutes. A sample was classified positive, when both the control and test line were visible within 10 minutes.

Table 5.4. **Interference studies with the gabControl iFOBT test using clinically defined OB negative samples supplemented with various potential interfering substances.**

# Tests done	Clinically defined negative OB specimen supplemented with following substances	gabControl iFOBT Test result (visible test and control band on test cassette strip = positive) (no visible test band, but visible control band on test cassette strip = negative)
5	Acetamiophen, 20 mg/dl	All (-)
5	Acetyl salicylic Acid, 20 mg/dl	All (-)
5	Ascorbic Acid, 20 mg/dl	All (-)
5	Atropine, 20 mg/dl	All (-)
5	Bilirubin, 60 mg/dl	All (-)
5	Caffeine, 20 mg/dl	All (-)
5	Creatinine, 20 mg/dl	All (-)
5	Gentesic Acid, 20 mg/dl	All (-)
5	Glucose, 2000 mg/dl	All (-)

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In conclusion, none of the above tested substances showed any interferences with neither a clinically defined negative nor a positive specimen. Negative specimen samples with supplementation of potentially interfering substances gave consistently negative test results, whereas specimen samples positive to OB scored consistently positive.

#### 4.8 Conclusion:

The relative sensitivity of the gabControl iFOBT test is 96.8 % (60/62) and the relative specificity is 99.6% (685/688) when compared to the InSure iFOBT. It is important to note that InSure iFOBT is not the gold-standard for colorectal cancer (CRC).

Gabmed's iFOBT test has a sensitivity of 68% and a specificity of 99.1% relative to the colonoscopy (the gold-standard) for the diagnosis of colorectal cancer, significantly higher than the other leading iFOBT.

The absolute sensitivity of gabControl iFOBT test is 50ng/ml of hemoglobin, which is one of the highest sensitivity on the market today.

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16/25

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## 5 Read Time Studies

Read Time Study; Remove the test disk from the foil packet, and place it on a flat, dry surface. Holding the sample dropper above the test disk, squeeze 2 to 3 drops of specimen into the sample well. As the test begins to work, we see purple color move across the result window in the center of the test disk and calculate the time. The read time for control and test line is summarized as in Table 5-1.

**Table 5-1: Read Time Study**

	Specimen	Control Line (in minutes, average)	Test Line (in minutes, average)
5	negative	+ (3)	N
5	Weak positive	+ (3)	P (10)
5	positive	+ (3)	P (5)
5	Strong Positive	+ (3)	P (2)

Conclusion; Read time for Control Line are 3 minutes and for Test Line are 5 to 10 minutes.

The Read Time Studies were done at laboratory of our company.

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17/25

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## 6 Reproducibility Study

The reproducibility study was carried in three different sites. Negative, weak positive and positive specimen samples had been tested by total 120 gabControl iFOBT test kit. The samples were tested two times in the same day, and in two different assays, each day for 20 days. The study results are shown in Table 6-1.

**Table 6-1: Reproducibility Study**

Days	OB (Hemoglobin)		
	Negative	Weak Positive (50ng/ml)	Strong Positive (200ng/ml)
	<b>Test Result (number of tests)</b>		
1	N (2)	P (2)	P (2)
2	N (2)	P (2)	P (2)
3	N (2)	P (2)	P (2)
4	N (2)	P (2)	P (2)
5	N (2)	P (2)	P (2)
6	N (2)	P (2)	P (2)
7	N (2)	P (2)	P (2)
8	N (2)	P (2)	P (2)
9	N (2)	P (2)	P (2)
10	N (2)	P (2)	P (2)
11	N (2)	P (2)	P (2)
12	N (2)	P (2)	P (2)
13	N (2)	P (2)	P (2)
14	N (2)	P (2)	P (2)
15	N (2)	P (2)	P (2)
16	N (2)	P (2)	P (2)
17	N (2)	P (2)	P (2)
18	N (2)	P (2)	P (2)
19	N (2)	P (2)	P (2)
20	N (2)	P (2)	P (2)

Conclusion: The kits were tested two times in the same day, and in two different assays each day for 20 days. This permits separate tests of between day, between-assay and within-day which they are shown consistent results.

The reproducibility study was performed at three different sites. The test results are shown Table 6-2.

**Table 6-2: Reproducibility Study Based on Study Sites**

Test sites	OB (Hemoglobin)		
	Negative	Weak Positive (50ng/ml)	Strong Positive (200ng/ml)
	<b>Test Readings ( number of tests)</b>		
1	N (40)	P (40)	P (40)
2	N (40)	P (40)	P (40)
3	N (40)	P (40)	P (40)

Conclusion: The tests results of three sites were in 100% agreement.

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### 7.1 All Materials;

#### 7.1.1 Materials;

Absorbent Filter Paper Type I  
 Specifications: Wet strength: 70  
 Composition: 100% cotton fiber  
 Absorbent Filter Paper Type II  
 Specifications: Wicking time: 3 minutes  
 Composition: 100 % glass fiber

#### 7.1.2 PVC;

Specification: membrane: cellulose nitrate  
 Backing: polyester film  
 Pore size: 5 um to 10 um  
 Thickness: 125 um  
 Wicking rate: < 30 sec/20 mm

#### 7.1.3 Major Antibodies;

Test zone: Mouse anti-hemoglobin antibody, 6 ug/test  
 Anti-hemoglobin antibody-Colloidal Gold Conjugate, 15 ug/test  
 Control zone: Goat anti-(Mouse IgG) Polyclonal antibody, 4 ug/test  
 Sample Collection Device Plastic tube 15 grams/ 10 cc

7.1.4 Buffer Solution Size 2 ml

#### 7.1.5 Major Buffers;

PBS: NaCl 0.1%

### 7.2 Production Process:

a. Nitrocellulose membrane is coated with a Goat Anti – (Mouse IgG) Polyclonal Antibody solution in the Control Reaction Zone, and coated with an anti-OB antibody solution in the Test Reaction Zone. The coated membrane is dried overnight, blocked with a PBS/Triton X-100 buffer for 5 minutes, and then dried over night again.

b. A buffered solution of anti-OB antibody – Colloidal Gold Conjugate is sprayed on fiberglass and then lyophilized in a 48-hour cycle. The fiberglass is trimmed to within 5 mm of the lyophilized conjugate band.

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19/25

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c. The blocked and dried membrane is then applied to an adhesive-backed vinyl strip. An absorbent material is applied to the vinyl strip so that it is in contact with the membrane close to the Control Reaction Zone. The lyophilized conjugate – coated fiberglass is applied to the vinyl strip are then die-cut at a 4 mm width.

d). The die-cut strips are assembled in plastic housings and sealed in moisture-proof foil pouches along with a desiccant packet and a disposable plastic pipette.

### **7.3 LOT # assignment**

Lot number is assigned when the QC is passed for the LOT in according to absolute sensitivity of 50ng/ml of hemoglobin.

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20/25

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## 7.4 Flow chart of the manufacturing process


Figure 7-1: Manufacturing Process Flow Chart

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## 7.5 Results of Stability studies

### 7.5.1 The Accelerating Stability Study;

The test also should pass following accelerating stability study in order assigning an Exp. date. The stability of test is depending on the temperature of storage room. The antibodies of the test kit will die if keeping them at high temperature room. The test kits are no longer working after antibodies dead. Our experience of the relationship between shelf life time and room temperature is presenting at following tab.

**Table 7-1: Accelerated Stability Study**

DESIGNED PRODUCT STORAGE: 4-30 C	
Months	Temperature
	56 C
6	1 Weeks
12	2 Weeks
18	3 Weeks
30	5 Weeks

We kept 25 tests in an oven at 56 C for 3 weeks. The accelerated test kit compared with normal test kit testing for negative and positive control. They show identical response in term of response time and results of the test. The conclusion: 5 weeks acceleration test indicated that they have an at least 30 month shelf life. The expiration was selected as 30 months from date of manufacturing.

Each LOT samples are also retained for real time stability test. The results indicated the test kits are normal function within expiration day and are consistent with the accelerated study.

### 7.5.2 Real Time Stability Study;

#### 7.5.2.1 Method

Reproducibility and stability studies were carried out on three different lots of the gabControl iFOBt Test. The lots, upon approval (t=0) were stored at 15 – 38 degrees C, and tested over an eighteen month time period from date of manufacture. In-house controls were tested according to the Instructions for use with the three lots at regular intervals (0, 3,6,9,12,15, 18 and 30 months from date of approval).

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## 7.5.2.2 Acceptance Criteria

Each stability test is passed if:

- The test gave a negative result when tested with a OB-negative specimen sample.
- The test gave a positive result when tested with a OB-positive (hemoglobin = 50ng/ml) specimen sample.

Results are shown in the following table.

**Table 7-2: Real Time Stability Study (3 Lots)**

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The test results show that all three LOTS performed consistently during 30 months of real time stability study.

Verteiler Original

Produktordner

Revision/ Fassung Nr.: 001

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23/25

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## 8 References

1. Walsh JM, Terdiman JP. Colorectal cancer screening: scientific review [review]. *Jama* 2003;289:1288-96.
2. Rex DK. Current colorectal cancer screening strategies: overview and obstacles to implementation [review]. *Rev Gastroenterol Disord* 2002;2 Suppl 1:S2-11.
3. Ahlquist DA, Shuber AP. Stool screening for colorectal cancer: evolution from occult blood to molecular markers. *Clin Chim Acta* 2002;315:157-68.
4. Otto S, Eckhardt S. Early detection for colorectal cancer: new aspects in fecal occult blood screening. *J Surg Oncol* 2000;75:220-6.
5. Mandel JS, Bond JH, Church TR, Snover DC, Bradley GM, Schuman LM, et al. Reducing mortality from colorectal cancer by screening for fecal occult blood. Minnesota Colon Cancer Control Study. *N Engl J Med* 1993;328:1365-71.
6. Mandel JS, Church TR, Ederer F, Bond JH. Colorectal cancer mortality: effectiveness of biennial screening for fecal occult blood. *J Natl Cancer Inst* 1999;91:434-7.
7. Hardcastle JD, Chamberlain JO, Robinson MH, Moss SM, Amar SS, Balfour TW, et al. Randomized controlled trial of fecal-occult-blood screening for colorectal cancer. *Lancet* 1996;348:1472-7.
8. Scholefield JH, Moss S, Sufi F, Mangham CM, Hardcastle JD. Effect of faecal occult blood screening on mortality from colorectal cancer: results from a randomized controlled trial. *Gut* 2002;50:840-4.
9. Kronborg O, Fenger C, Olsen J, Jorgensen OD, Sondergaard O. Randomised study of screening for colorectal cancer with fecal-occult-blood test. *Lancet* 1996;348:1467-71.
10. Jorgensen OD, Kronborg O, Fenger C. A randomised study of screening for colorectal cancer using faecal occult blood testing: results after 13 years and seven biennial screening rounds. *Gut* 2002;50:29-32.
11. Rasmussen M, Fenger C, Kronborg O. Diagnostic yield in a biennial Hemoccult-II screening program compared to a once-only screening with flexible sigmoidoscopy and Hemoccult-II. *Scand J Gastroenterol* 2003;38:114-8.

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24/25

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## 9 Version History

Product	changes/ revision	reason	released
FOB M04K11	1.0		April 2008

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actim™

## Actim Fecal Blood evaluations

### Actim Fecal Blood comparison with Hemocult test

Presented at: Pediatric Congress with surgery problematics, Czech Republic by Kocna P. M.D and PhD Dvorák M. M.D. Ústav klinické biochemie a 4. interí klinika VFN a1.LF UK Praha.

There were 213 completely evaluated patients. Totally 280 patients were drafted to the study. Colonoscopy was performed with totally 218 patients. Actim Fecal Blood test was not performed with 5 patients.

Colorectal disease	Actim Fecal Blood		Hemocult	
	+	-	+	-
+	49	40	23	66
-	8	116	5	119

	Actim Fecal Blood	Hemocult
<b>Specificity (%)</b>	93.55	95.97
<b>Sensitivity (%)</b>	55.01	25.84
<b>Precision (%)</b>	77.46	66.67
<b>PPV (%)</b>	85.96	82.14
<b>NPV (%)</b>	74.36	64.32

The sensitivity of the Actim Fecal Blood is clearly better than the sensitivity of the Hemocult test. However, the sensitivity of neither of the two tests is very good. There were only a few false positives and therefore the specificity of both tests is good.

actim™

## Actim Fecal Blood evaluations

### Evaluation at Tampere University Hospital, Finland

Two consecutive fecal samples were collected from 41 patients and were analyzed by three alternative test methods (Fecatwin, Hemolex and Actim Fecal Blood). Gastroendoscopy and colonoscopy were then performed. The findings were classified as colorectal disease or disease in the upper part of the gastrointestinal tract or negative. Three patients were simultaneously diagnosed with upper gastrointestinal tract and colorectal disease.

The sensitivity and specificity of each test was calculated for Colorectal disease (including the three patients with also a disease of the upper gastrointestinal tract) and for Upper bleeding separately.

Colorectal disease	Actim Fecal Blood		Hemolex		Fecatwin	
	+	-	+	-	+	-
+	8	3	6	5	7	4
-	0	12	0	12	4	8
Upper Bleeding						
+	5	8	3	10	4	9
-	0	12	0	12	4	8

Colorectal disease	Actim Fecal Blood	Hemolex	Fecatwin
<b>Sensitivity</b>	72.7%	54.5 %	63.6 %
<b>Specificity</b>	100 %	100 %	66.7 %
Upper Bleeding			
<b>Sensitivity</b>	38.5 %	23.1 %	30.8 %
<b>Specificity</b>	100 %	100 %	66.7 %

Note: The Colorectal disease group includes the three cases with simultaneous diagnosis of upper bleeding and colorectal disease

The Actim Fecal Blood test had the highest sensitivity in both groups. No false positive results were detected for patients without pathological findings in the endoscopies. The specificity of the two immunological tests (Hemolex and Actim Fecal Blood) was 100%, but with the Guaiac test (Fecatwin) it was only 66,7%.

Actim Fecal Blood test had the best combination of sensitivity and specificity for detecting colorectal disease. For bleeding from the upper GI tract, the specificity of the immunological tests is good and there were no false positive results. However, there were quite a many false negative results and therefore the sensitivity is low.

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